

STEREO ACCESSORIES



outer connected illumination
12V/10W



Light source 24W/150W



Square illumination
110V/9W/220/9W



RFL10T



RFL8



RFL-10-1



RFL8-1



RFL10



LED-48T



LED-48



LED-60T



LED-96T

Model Parameters	Single Light Guide	Dual Light Guide	Ring Light Guide	
	Length (mm)	550	550	900
Optical Fiber (mm)	Φ5.6	Φ5.6	Φ10	Φ8
Ring Inner Dia. (mm)	/	/	Φ55.5	Φ55.5

Fluorescent ring lamp:

Model	RFL8	RFL8-1	RFL8-2	RFL10-1	RFL10-2	RFL10-2T
Voltage	110V / 220V					
Power	8W			10W		
Color Temperature	2700K / 6000K					
Light Adjustment	No light adjustment					Light adjustment

LED ring lamp

Model	Connection mode between light source and light- regulating part	Be allocated with ballast (input 110V/220V)	Luminous fLux density (working distance 100mm)
LED-24	No light adjustment	4.5V / 1A	4500
LED-48			9000
LED-24T	Integrated light adjustment	24V / 250mA	520-2400

LED-48T			500-18000			
LED-60T	Detachable light adjustment	12V/ 600mA	Four variable lighting modes			
			Fully bright	1/4	1/2	3/4
LED-96T		12V	6500			
			7step regulation 2600- 25000	8000	15000	21000



Selective accessories	Technical parameters	Selective accessories	Technical parameters
1X Objective	W.D.57mm	WF5X Eyepiece	Φ20mm
2X Objective	W.D.80mm	WF10X Eyepiece	Φ20mm
3X Objective	W.D.61mm	WF15X Eyepiece	Φ15mm
4X Objective	W.D.57mm	SWF10X Eyepiece	Φ22mm
5X Objective	W.D.66mm	SWF20X Eyepiece	Φ10mm



TD- II Universal Stand with Large Base



TD- I Universal Stand with Large Base



Darkfield attachment



(I) Jewel tweezers (II)